

Aristotle University of Thessaloniki



Physics Department

Exercise 1

Spectrum analysis on ^{55}Fe using ROOT

Author: Konstantinos E. Stergiou

April 15, 2025

Contents

Complete ROOT Script	3
Step 1: Histogram creation	5
Step 2: Data loading	5
Step 3: Fit function definition	5
Step 4: Initial parameters	5
Step 5: Fit	5
Step 6: Plot	6
Step 7: Peak ratio	6
Results	7

Introduction

In the exercise below, we analyze a gas detector spectrum obtained by irradiation with 5.9 keV X-rays from a ^{55}Fe source. Due to electronic noise, the escape peak at ~ 2.7 keV is partially masked. A multi-Gaussian fit with linear background is performed to extract and quantify the relative peak intensities.

Key points

- **Source:** ^{55}Fe emits characteristic X-rays at 5.9 keV (85%) and 6.5 keV (15%)
- **Escape Peak:** Argon K-shell escape reduces energy by ~ 3.2 keV, forming a peak at ~ 2.7 keV
- **Method:** Fit includes 3 Gaussians + linear background
- **Goal:** Determine the intensity ratio between the 5.9 keV and escape peaks

Complete ROOT Script

Listing 1: Complete ROOT script for ^{55}Fe spectrum analysis

```
1 // analyze_55Fe.C
2 // ROOT script for 55Fe spectrum analysis
3
4 void analyze_55Fe() {
5     // 1. Create histogram
6     TH1F *spc4 = new TH1F("spc4", "55Fe Spectrum;Channel;
7         Counts", 1024, 0, 1024);
8
9     // 2. Read data file
10    FILE *fp4 = fopen("QMO_5310.TXT", "r");
11    if (!fp4) {
12        std::cerr << "Error: Cannot open file!" << std::endl
13        ;
14        return;
15    }
16
17    char line[80];
18    int chnl, entr, nn;
19    while (fgets(line, 80, fp4)) {
20        sscanf(line, "%d,%d %d", &chnl, &entr, &nn);
21        spc4->SetBinContent(chnl, entr);
22    }
23    fclose(fp4);
24
25    // 3. Define fit function (3 Gaussians + linear
26    background)
27    TF1 *fitFunc = new TF1("fitFunc",
28        "[0]*exp(-0.5*((x-[1])/[2])^2) + " // 5.9 keV peak
29        "[3]*exp(-0.5*((x-[4])/[5])^2) + " // 6.5 keV peak
30        "[6]*exp(-0.5*((x-[7])/[8])^2) + " // Escape peak
31        (~2.7 keV)
32        "[9] + [10]*x", // Linear
33        background
34        0, 1024);
35
36    // 4. Initial parameters (adjust based on your data)
37    fitFunc->SetParameters(
38        1120, 220, 30, // 5.9 keV peak (amplitude, mean,
39        sigma)
40        350, 430, 100, // 6.5 keV peak
41        5850, 180, 20, // Escape peak
42        -3, 0.1 // Background (constant + slope)
43    );
44
45    // 5. Perform the fit
46    spc4->Fit("fitFunc", "R"); // "R" for fit range
47
48    // 6. Plot results
```

```

43   TCanvas *c2 = new TCanvas("c2", "55Fe Fit Results",
44       1000, 800);
45   spc4->Draw();
46   fitFunc->Draw("same");
47
48   // Draw individual components
49   TF1 *g1 = new TF1("g1", "gaus", 0, 1024);
50   g1->SetParameters(fitFunc->GetParameter(0), fitFunc->
51       GetParameter(1), fitFunc->GetParameter(2));
52   g1->SetLineColor(kRed);
53   g1->Draw("same");
54
55   TF1 *g2 = new TF1("g2", "gaus", 0, 1024);
56   g2->SetParameters(fitFunc->GetParameter(3), fitFunc->
57       GetParameter(4), fitFunc->GetParameter(5));
58   g2->SetLineColor(kBlue);
59   g2->Draw("same");
60
61   TF1 *g3 = new TF1("g3", "gaus", 0, 1024);
62   g3->SetParameters(fitFunc->GetParameter(6), fitFunc->
63       GetParameter(7), fitFunc->GetParameter(8));
64   g3->SetLineColor(kGreen);
65   g3->Draw("same");
66
67   // 7. Calculate peak ratios
68   double area_59keV = g1->Integral(g1->GetParameter(1)-3*
69       g1->GetParameter(2),
70       g1->GetParameter(1)+3*g1->
71       GetParameter(2));
72   double area_escape = g3->Integral(g3->GetParameter(1)-3*
73       g3->GetParameter(2),
74       g3->GetParameter(1)+3*g3->
75       GetParameter(2));
76
77   std::cout << "\nResults:\n";
78   std::cout << "5.9 keV peak area: " << area_59keV << "\n"
79       ;
80   std::cout << "Escape peak area: " << area_escape << "\n"
81       ;
82   std::cout << "Ratio (5.9keV/Escape): " << area_59keV/
83       area_escape << "\n";
84
85   // 8. Save results
86   c2->SaveAs("55Fe_fit_results.png");
87 }

```

Step 1: Histogram creation

```
1 // 1. Create histogram
2 TH1F *spc4 = new TH1F("spc4", "55Fe Spectrum;Channel;Counts"
    , 1024, 0, 1024);
```

Explanation: Creates a 1D histogram (TH1F) with:

- 1024 bins (channels) from 0 to 1024
- Title "55Fe Spectrum"; x-axis "Channel"; y-axis "Counts"
- spc4 is the pointer to the histogram object

Step 2: Data loading

```
1 // 2. Read data file
2 FILE *fp4 = fopen("QMO_5310.TXT", "r");
3 // Read and parse lines, fill histogram
```

Data from the detector is read from a '.TXT' file line-by-line and each entry is loaded into the histogram.

Step 3: Fit function definition

```
1 // 3. Define fit function (3 Gaussians + linear background)
2 TF1 *fitFunc = new TF1("fitFunc",
3     "[0]*exp(-0.5*((x-[1])/[2])^2) + "
4     "[3]*exp(-0.5*((x-[4])/[5])^2) + "
5     "[6]*exp(-0.5*((x-[7])/[8])^2) + "
6     "[9] + [10]*x",
7     0, 1024);
```

The fit function consists of 3 Gaussians and a linear background model.

Step 4: Initial parameters

```
1 // 4. Initial parameters
2 fitFunc->SetParameters(
3     1120, 220, 30,      // 5.9 keV
4     350, 430, 100,     // 6.5 keV
5     5850, 180, 20,     // Escape
6     -3, 0.1            // Background
7 );
```

Estimated initial values for peak amplitudes, means and widths, based on expected peak locations.

Step 5: Fit

```
1 // 5. Perform the fit
2 spc4->Fit("fitFunc", "R");
```

Performs the fit using the specified function over the full histogram range.

Step 6: Plot

```
1 // 6. Plot results
2 TCanvas *c2 = new TCanvas("c2", "55Fe Fit Results", 1000,
   800);
3 spc4->Draw();
4 fitFunc->Draw("same");
```

The total fit and each component Gaussian are drawn on a ROOT canvas.

Step 7: Peak ratio

```
1 // 7. Calculate peak ratios
2 double area_59keV = g1->Integral(...);
3 double area_escape = g3->Integral(...);
4 double ratio = area_59keV / area_escape;
```

Using numerical integration over $\mu \pm 3\sigma$, the area under each peak is computed and their ratio is calculated.

Results




Figure 1: Fit results for the ^{55}Fe spectrum showing the three Gaussian components (red: 5.9 keV, blue: 6.5 keV, green: escape peak) and linear background.

The peak areas computed from the fitted Gaussians are:

- **5.9 keV peak area:** 87,295.7
- **Escape peak area:** 344,669

The ratio of the main 5.9 keV peak to the escape peak is:

$$\text{Ratio}_{5.9 \text{ keV/Escape}} = \boxed{0.253}$$

This ratio characterizes the detector's response to the 5.9 keV X-rays versus the argon escape events, as visualized in Fig. 1.